Gebrauchsanleitung | Operating manual | Mode  
d'emploi | Instrucciones de manejo | Istruzione | 操作  
手册

# Transferpette® electronic

Mikroliterpipetten | Micropipettes

# Impressum

BRAND GMBH + CO KG

Otto-Schott-Str. 25

97877 Wertheim (Germany)

T +49 9342 808 0

F +49 9342 808 98000

info@brand.de

www.brand.de

## Languages

Gebrauchsanleitung.....	2
Operating manual.....	45
Mode d'emploi .....	82
Instrucciones de manejo .....	120
Istruzione.....	158
操作手册 .....	197

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## EU-Conformity Declaration

Handelsname des Produktes / Trade name of the device:	Kolbenhubpipette Transferpette®/electronic, Piston-operated pipette Transferpette® electronic
Produktvarianten / Device variations:	-8, /-12, (0,5-10 µl, 1-20 µl, 2-20 µl, 5-100 µl, 10-200 µl, 15-300 µl, 50-1.000 µl, 50-1250 µl, 250-5.000 µl) -8, /-12, (0,5-10 µl, 1-20 µl, 2-20 µl, 5-100 µl, 10-200 µl, 15-300 µl, 50-1.000 µl, 50-1250 µl, 250-5.000 µl)
Art.-Nr. / Cat.-No.:	705299, 705300, 705303, 705306, 705307, 705309, 705310, 705313, 705316, 705317, 705319, 705320, 705323, 705326, 705327, 705329, 705330, 705333, 705336, 705337, 705339, 705340, 705343, 705346, 705347, 705398, 705399, 705400, 705403, 705404, 705406, 705409, 705410, 705413, 705414, 705416, 705419, 705420, 705423, 705424, 705426, 705429, 705430, 705433, 705434, 705436, 705448, 705449, 705450, 705453, 705454, 705456, 705459, 705460, 705463, 705464, 705466, 705469, 705470, 705473, 705474, 705476, 705479, 705480, 705483, 705484, 705486
Hersteller / Manufacturer:	BRAND GMBH + CO KG Otto-Schott-Str. 25 97877 Wertheim

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# UK Declaration of Conformity

Trade name of the device: Piston-operated pipette Transferpette<sup>®</sup>-electronic

Device variations: -8, /-12, (0,5-10 µl, 1-20 µl, 2-20 µl, 5-100 µl, 10-200 µl, 15-300 µl, 50-1.000 µl, 50-1250 µl, 250-5.000 µl)

Cat.-No.: 705299, 705300, 705303, 705306, 705307, 705309, 705310, 705313, 705316, 705317, 705319, 705320, 705323, 705326, 705327, 705329, 705330, 705333, 705336, 705337, 705339, 705340, 705343, 705346, 705347, 705398, 705399, 705400, 705403, 705404, 705406, 705409, 705410, 705413, 705414, 705416, 705419, 705420, 705423, 705424, 705426, 705429, 705430, 705433, 705434, 705436, 705448, 705449, 705450, 705453, 705454, 705456, 705459, 705460, 705463, 705464, 705466, 705469, 705470, 705473, 705474, 705476, 705479, 705480, 705483, 705484, 705486

Manufacturer: BRAND GMBH + CO KG  
Otto-Schott-Str. 25  
97877 Wertheim · Germany

We, as the manufacturer of the above described device(s) take sole responsibility for and hereby declare that the described device(s) meet the provisions of the following Regulation(s)/Directive(s):	Applied harmonized standards:
The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012	EN IEC 63000:2018
Electromagnetic Compatibility Regulations 2016 (SI 2016 No. 1091, as amended by SI 2019 No. 696)	EN 61326-1:2013
<b>Other applied standards:</b>	<b>Scope:</b>
EN 61010-1:2010, EN 61010-1:2010/A1:2019/AC:2019-04, EN 61010-1:2010/A1:2019	Electrical and functional safety

Wertheim December 15, 2021.

#08.19.01.01

*P. Ziemeck*

Patrick Ziemeck  
Managing Director Technology

*i.A. Siegfried Ott*

i.A. Siegfried Ott  
Regulatory Affairs

This document declares the accordance with the named regulations, but does NOT assure specific properties

BRAND GMBH + CO KG - info@brand.de- www.brand.de

# Declaration of conformity – China RoHS 2

BRAND GMBH + CO KG has made reasonable efforts to ensure that hazardous materials and substances may not be used in BRAND products.

In order to determine the concentration of hazardous substances in all homogeneous materials of the subassemblies, a "Product Conformity Assessment" (PCA) procedure was performed. As defined in GB/T 26572 the "Maximum Concentration Value" limits (MCV) apply to these restricted substances:

+ Lead (Pb):	0.1%	+ Hexavalent chromium (Cr(+VI)):	0.1%
+ Mercury (Hg):	0.1%	+ Polybrominated biphenyls (PBB):	0.1%
+ Cadmium (Cd):	0.01%	+ Polybrominated diphenyl ether (PBDE):	0.1%

## Environmental Friendly Use Period (EFUP)

EFUP defines the period in years during which the hazardous substances contained in electrical and electronic products will not leak or mutate under normal operating conditions. During normal use by the user such electrical and electronic products will not result in serious environmental pollution, cause serious bodily injury or damage to the user's assets.

The Environmental Friendly Use Period for BRAND instruments is 40 years.

此表格是按照SJ/T 11364-2014中规定所制定的。  
This table is created according to SJ/T 11364-2014.

## Material content declaration for BRAND products

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**Note:** Table applies to all products. Some of the components or parts listed above may not be part of the enclosed product.

- O: 表示该有毒有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下。
- O: Indicates that the above mentioned hazardous substance contained in all homogeneous materials of the part is below the required limit as defined in GB/T 26572.
- X: 表示该有毒有害物质至少在该部件某一均质材料中的含量超出GB/T 26572规定的限量要求。
- X: Indicates that the above mentioned hazardous substance contained in at least one of the homogeneous materials of this part is above the required limit as defined in GB/T 26572.

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Apart from the disclosures in the above table, the subassemblies are not intentionally manufactured or formulated with lead (Pb), mercury (Hg), cadmium (Cd), hexavalent chromium (Cr+VI), polybrominated biphenyls (PBB), and polybrominated diphenyl ethers (PBDE).

Products manufactured by BRAND may enter into further devices or can be used together with other appliances. With these third party products and appliances in particular, please note the EFUP labeled on these products. BRAND will not take responsibility for the EFUP of those products and appliances.

Place, date: Wertheim, 17|05|2021

**Patrick Ziemeck**  
(Managing Director  
Technology)

**BRAND GMBH + CO KG** | Otto-Schott-Str. 25 | 97877 Wertheim | Germany  
T +49 9342 808-0 | F +49 9342 808-98000 | info@brand.de | www.brand.de

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#### Indien:

BRAND Scientific Equipment Pvt. Ltd.  
303, 3rd Floor, 'C' Wing, Delphi

#### China:

BRAND (Shanghai) Trading Co., Ltd.  
Guangqi Culture Plaza  
Room 506, Building B

Hiranandani Business Park,  
Powai  
Mumbai-400 076 (India)  
T +91 22 42957790  
F +91 22 42957791  
info@brand.co.in  
www.brand.co.in

No. 2899, Xietu Road  
Shanghai 200030 (P.R. China)  
T +86 21 6422 2318  
F +86 21 6422 2268  
info@brand.com.cn  
www.brand.cn.com

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
# Table of contents

<b>1 Introduction</b> .....	<b>46</b>	10.4 Rechargeable battery regeneration function.....	72
1.1 Scope of delivery.....	46		
1.2 Terms of use.....	46		
<b>2 Safety regulations</b> .....	<b>47</b>	<b>11 Troubleshooting</b> .....	<b>73</b>
2.1 General safety regulations.....	47	<b>12 Product markings</b> .....	<b>74</b>
2.2 Intended use.....	47	<b>13 Order Information</b> .....	<b>75</b>
2.3 Limitations of use.....	48	13.1 Ordering Information .....	75
2.4 Application restrictions.....	48	13.2 Spares .....	77
2.5 Operating exclusions .....	48	13.3 Additional accessories .....	78
2.6 Rechargeable battery and power adapter specifications .....	48	13.4 Accessories for Transferpette® electronic 10–200 µl and 50–1000 µl .....	78
<b>3 Functional and operational components</b> <b>49</b>		<b>14 Repairs</b> .....	<b>79</b>
<b>4 Commissioning</b> .....	<b>50</b>	14.1 Sending for repair.....	79
4.1 First Steps.....	50	<b>15 Calibration service</b> .....	<b>80</b>
4.2 Set volume.....	51	<b>16 Warranty</b> .....	<b>81</b>
4.3 Set the aspiration and dispensing speed .....	52	<b>17 Disposal</b> .....	<b>81</b>
4.4 Pipetting.....	53		
4.5 Directly trigger blow-out.....	54		
<b>5 Pipetting programs</b> .....	<b>55</b>		
5.1 PIP mode .....	55		
5.2 PIPmix mode .....	57		
5.3 revPIP mode .....	58		
5.4 Electrophoresis (GEL) mode .....	60		
5.5 DISP mode .....	61		
<b>6 Checking the volume</b> .....	<b>64</b>		
<b>7 Accuracy table</b> .....	<b>65</b>		
<b>8 Adjustment – Easy Calibration</b> .....	<b>66</b>		
8.1 Adjustment .....	66		
8.2 Restore factory setting.....	67		
<b>9 Disinfection/autoclaving</b> .....	<b>68</b>		
9.1 UV sterilization .....	68		
9.2 Autoclaving.....	68		
9.3 Reference run (rEF) .....	68		
<b>10 Maintenance</b> .....	<b>69</b>		
10.1 Disassembling/cleaning (up to 1,000 µl) .....	69		
10.2 Disassembly/cleaning (250–5000 µl) .....	70		
10.3 Charge and replace rechargeable battery .....	71		

# 1 Introduction

## 1.1 Scope of delivery

Transferpette® electronic ,a rechargeable battery, the power adapter with battery charging cable, silicone grease, this operating manual, and one sample bag of pipette tips.

## 1.2 Terms of use

- Carefully read the operating manual before using the device for the first time.
- The operating manual is part of the device and must be kept in an easily accessible place.
- Be sure to include the operating manual if you transfer possession of this device to a third party.

### 1.2.1 Hazard levels

The following signal words identify possible hazards:

Signal word	Meaning
DANGER	Will lead to serious injury or death.
WARNING	May lead to serious injury or death.
CAUTION	May lead to minor or moderate injuries.
NOTICE	May lead to property damage.

### 1.2.2 Format

Format	Meaning	Format	Meaning
<b>1. Task</b>	Indicates a task.	>	Indicates a condition.
a., b., c.	Indicates the individual steps of a task.	⇒	Indicates a result.

### 1.2.3 Operating manual symbols

Icon	Meaning	Icon	Meaning
	Hazardous area		Explosion hazard

## 2 Safety regulations

### 2.1 General safety regulations

#### Please read carefully!

The instrument Transferpette® electronic can be used in combination with hazardous materials, work processes and equipment. However, the operating manual cannot cover all of the safety issues that may occur in doing so. It is the user's responsibility to ensure compliance with the safety and health regulations and to specify the corresponding restrictions before use.

1. Every user must read and understand this operating manual before operation.
2. Follow the general hazard instructions and safety regulations, e.g. wear protective clothing, eye protection and protective gloves.  
When working with infectious or hazardous samples, the standard laboratory rules and precautions must be adhered to.
3. Observe all specifications provided by the reagent manufacturer.
4. Do not operate the instrument in potentially explosive atmospheres and do not pipette highly flammable media.
5. Use the instrument only for pipetting liquids within the defined limitations and restrictions of use. Comply with the operating exclusions; see Operating exclusions, p. 48. If in doubt, contact the manufacturer or supplier.
6. Always perform work in a way that does not endanger the user or any other person. Avoid splashes. Use only suitable vessels.
7. Avoid touching the tip opening when working with aggressive media.
8. Never use force.
9. Use only original spare parts. Do not attempt to make any technical modifications.  
Do not disassemble the instrument any further than is described in the operating manual.
10. Always check that the instrument is in proper working condition before use. If there are any signs of the instrument malfunctioning (e.g., sluggish piston, leakage), stop pipetting immediately, and follow the instructions in the chapter Troubleshooting, p. 73. Contact the manufacturer, if necessary.
11. The original rechargeable battery may not be exchanged for batteries or rechargeable batteries from other manufacturers.
12. Use only the original power adapter to charge the nickel-metal hydride battery.
13. Protect the power adapter from moisture and use it only with this instrument.
14. Dispose of batteries according to the applicable regulations.

#### WARNING

#### Potential risk of explosion due to damaged battery

In extreme cases, improper handling of the instrument or the battery (short circuit, mechanical damage, overheating, etc.) may cause the battery to explode.

### 2.2 Intended use

Transferpette® electronic is a microprocessor-controlled, battery-operated piston-stroke pipette based on the air cushion principle for pipetting aqueous solutions of medium density and viscosity. If the instrument is handled correctly, the sample to be dispensed comes into contact only with the tip and not with the Transferpette® electronic.

## 2.3 Limitations of use

This instrument is intended for pipetting samples, within the following limitations:

- Operating temperature of instrument and reagent should be between +15 °C and +40 °C (59 °F to 104 °F) (other temperatures upon request)
- Vapor pressure up to 500 mbar
- Viscosity: 260 mPa s 260 mPa s

For viscous media, the speed must be adjusted if necessary.

## 2.4 Application restrictions

Viscous and wetting liquids may compromise volumetric accuracy. Volumetric accuracy may also be affected when pipetting liquids whose temperature deviates from the ambient temperature by more than  $\pm 1\text{ °C}/\pm 1.8\text{ °F}$ .

## 2.5 Operating exclusions

The user is responsible for checking the compatibility of the device with the intended application. The device cannot be used:

- for liquids that attack polypropylene or polyvinylidene fluoride, polycarbonate/polybutylene terephthalate, polyether ether ketone, FKM or EPDM (flexible replacement pipette shafts)

Avoid aggressive vapors (risk of corrosion)!

The hand grip is not autoclavable.

## 2.6 Rechargeable battery and power adapter specifications

### Rechargeable battery


Nickel metal hydride rechargeable battery with 3 cylindrical single cells size AAA, 3.6 V, 700 mAh

### Power adapter

Output voltage 6.5 V DC, 200 mA

## 3 Functional and operational components

- |           |   |           |                     |
|-----------|---|-----------|---------------------|
| <b>1</b>  | Charging socket                             | <b>2</b>  | Battery compartment |
| <b>3</b>  | Menu selection key                          | <b>4</b>  | Arrow key (+)       |
| <b>5</b>  | 'Enter' key Input confirmation/Power 'On'*) | <b>6</b>  | Arrow key (-)       |
| <b>7</b>  | Finger rest                                 | <b>8</b>  | Handle              |
| <b>9</b>  | Pipetting shaft                             | <b>10</b> | Tip cone            |
| <b>11</b> | Tip ejection key                            | <b>12</b> | Pipetting key       |
| <b>13</b> | Display                                     |           |                     |

\*) The instrument is switched on by pressing the 'Enter' key. By subsequently pressing the pipetting key, the instrument is ready for pipetting.

The Transferpette® electronic switches itself off automatically 10 minutes after the last operation (auto power off).

The Transferpette® electronic fits ergonomically in your hand. For absolutely effortless operation of the function keys, the hand position can be further optimized by means of the finger grip, which is height-adjustable via a screw.

# 4 Commissioning


## 4.1 First Steps

### 1. Insert rechargeable battery


- a. Open the cover of the battery compartment.
- b. Insert the rechargeable battery into the compartment. Make sure that the plug of the rechargeable battery is firmly inserted into the socket in the instrument.
- c. Replace the cover of the battery compartment and close it.

### 2. Activate instrument

The Transferpette® electronic automatically requests a reference run immediately after the battery is inserted. After pressing the pipetting key, the reference run is carried out, and the instrument is ready for pipetting.

The display shows the factory-set default pipetting mode (PIF) and the respective nominal volume.

The aspiration and dispensing speed are set to maximum at the factory.

The simple volume and speed setting is described on the following pages.

## 4.2 Set volume


The volume is factory-set to the respective nominal volume of the Transferpette® electronic and can be easily and quickly changed individually.


- a.** Pressing one of the arrow keys directly selects a volume.  
 ⇒ 'VOL' flashes.

- b.** Pressing the arrow key (-) decreases the volume. Pressing and holding the arrow key causes the volume to change quickly.  
 ⇒ 'VOL' continues to flash.

- c.** Pressing the arrow key (+) increases the volume. Pressing and holding the arrow key causes the volume to change quickly.  
 ⇒ 'VOL' continues to flash.

- d.** To confirm the volume selection, press the 'Enter' key.  
 ⇒ The display now shows the newly set volume. Here, for example, the display of the PIP mode set by default.


### NOTICE

By pressing the menu key, any setting procedure can be canceled. The display then jumps to the next setting option or back to the initial display.

## 4.3 Set the aspiration and dispensing speed


Aspiration and dispensing speed are separately adjustable. When the menu is called up, the last speed set in each case is displayed. Five speed levels are available in each case.

### Set aspiration speed


**a.** Pressing the menu key once briefly takes you to the aspiration speed menu.

⇒ 'Speed' flashes.

**b.** Press the arrow keys (+\-) to select the speed level (e.g. level 5).


⇒ 'Speed' continues to flash.

**c.** Press 'Enter' key.



⇒ The display returns to the basic state of the respective set mode. Here, for example the display of the standard PIP mode.

### Set dispensing speed

**a.** Pressing the menu key twice briefly takes you to the dispensing speed menu.

⇒ 'Speed' flashes.

- b.** Press the arrow keys (+\-) to select the speed level (e.g. level 2).


⇒ 'Speed' continues to flash.

- c.** Press 'Enter' key.

⇒ The display returns to the basic state of the respective set mode. Here, for example the display of the standard PIP mode.

## 4.4 Pipetting


The volume is factory-set to the respective nominal volume of the Transferpette® electronic and can be easily and quickly changed individually, see Set volume, p. 51.


- a.** Attach tip vertically:  
Use the correct tips, in accordance with the volume range or color code!  
Make sure that the tips are firmly in place and leak tight. When using the flexible pipette shaft, attach an alternative ejector adjustment clip if necessary.  
Pipette tips are disposable products!
- b.** Aspirate liquid: Hold the instrument vertically, and immerse the tip 2–3 mm in the liquid. The liquid is aspirated by pressing the pipetting key. The arrow in the display points upwards (aspiration).

Leave the tip immersed in the liquid for a few seconds, so that the set volume is aspirated completely. This is especially important when pipetting viscous media and when using pipettes with large volumes.

Volume range	Immersion depth	Wait time
0.5 – 100 µl	2 – 3 mm	1 s
100 – 300 µl	2 – 4 mm	1 s
> 1000 µl	3 – 6 mm	3 s

- c. Dispense liquid: When the liquid aspiration is finished, the arrow in the display points downwards (dispensing). Place the pipette tip against the vessel wall. Hold the pipette at an angle of 30–45°. By pressing the pipetting key again, the liquid is dispensed completely with automatic blow-out. While doing this, wipe the pipette tip against the vessel wall.


- d. Ejecting a tip: Hold the pipette shaft over a suitable disposal container and press the tip ejection key.

### NOTICE


ISO 8655 requires that pipette tips are pre-wetted once before the actual pipetting procedure.

## 4.5 Directly trigger blow-out

The blow-out can also be triggered directly at any time if necessary.

- a. Call up the blow-out function: Press 'Enter' key. The display shows 'blo' for blow-out.

- b. Trigger blow-out: Press the pipetting key once to trigger the blow-out. The display will then jump back to the set pipetting mode (start position).

### NOTICE

During blow-out, the piston moves completely downwards. Ensure that any residual liquid is discharged safely. Keeping the pipetting key pressed keeps the piston down and thus prevents accidental aspiration of liquid. Releasing causes the piston to return to the start position.

## 5 Pipetting programs

Activity	Designation	Info
Normal pipetting	PIPmode, see PIP mode, p. 55	Standard program. A previously entered volume is aspirated and dispensed again
Pipetting for electrophoresis	GELmode, see Electrophoresis (GEL) mode, p. 60	Program for loading electrophoresis gels. A previously defined sample volume is aspirated at a high, variable speed and slowly released again.
Mixing samples	PIPMix mode, see PIPmix mode, p. 57	Program for mixing liquids. Sample is aspirated and dispensed repeatedly.
Reverse pipetting	revPIP mode, see revPIP mode, p. 58	Program especially for pipetting liquids with high viscosity, high vapor pressure, or foaming media.
Dispensing	DISPmode, see DISP mode, p. 61	Program for dispensing liquids. A volume taken up is dispensed again in partial steps.

### NOTICE


#### GEL mode

The GEL mode is not available for Transferpette® electronic 1000 µl and 5000 µl.


### 5.1 PIP mode

The default program – a previously entered volume is aspirated and dispensed again.


Volume and speed setting, see Set volume, p. 51 and Set the aspiration and dispensing speed, p. 52.


- a.** Call up menu selection: Pressing the menu key three times takes you to the program selection.  
⇒ 'Mode' flashes.


- b.** Set PIP mode: Use one of the arrow keys to scroll through the modes until 'PIP' appears.  
⇒ 'Mode' continues to flash.


- c.** Confirm PIP mode: Press 'Enter' key.  
 ⇨ The display now shows 'blo' for blow-out.

- d.** Preparing for pipetting: Press the pipetting key once to move the piston to its start position.  
 ⇨ The arrow in the display points upwards (aspiration).

- e.** Aspirate liquid: Press the pipetting key once to aspirate the liquid.

- f.** Dispense liquid: The liquid is dispensed by pressing the pipetting key once.  
 ⇨ The arrow in the display points downwards (dispensing).


- g.** Trigger blow-out? You do not have to do anything. When pipetting in PIP mode, the blow-out is performed automatically.


## 5.2 PIPmix mode

Program for mixing liquids. Sample is aspirated and dispensed repeatedly.


Volume and speed setting, see Set volume, p. 51 and Set the aspiration and dispensing speed, p. 52.


- a.** Call up menu selection: Pressing the menu key three times takes you to the program selection.  
 ⇒ 'Mode' flashes.


- b.** Set PIPmix mode: Use one of the arrow keys to scroll through the modes until 'PIPmix' appears.  
 ⇒ 'Mode' continues to flash.


- c.** Confirm PIPmix mode: Press 'Enter' key.  
 ⇒ The display now shows 'blo' for blow-out.


- d.** Preparing for pipetting: Press the pipetting key once to move the piston to its start position.  
 ⇒ The arrow in the display points upwards (aspiration).

- e.** Aspirate liquid: Press the pipetting key once to aspirate the liquid.

f. Dispense liquid in PIPmix mode: Pressing and holding the pipetting key causes the liquid to be alternately dispensed and aspirated. The display alternately shows the arrow symbol for aspiration or dispensing as well as the number of cycles.


g. End pipetting: Pressing the pipetting key once dispenses the liquid and triggers the blow-out. After dispensing the residual liquid (blow-out), the display jumps back to the set mode (start position).

**NOTICE**


The display shows a maximum of 19 cycles.

### 5.3 revPIP mode


Program especially for pipetting liquids with high viscosity, high vapor pressure, or foaming media. Volume and speed setting, see Set volume, p. 51 and Set the aspiration and dispensing speed, p. 52.







a. Call up menu selection: Pressing the menu key three times takes you to the program selection.  
 ⇨ 'Mode' flashes.

b. set revPIP mode: Use one of the arrow keys to scroll through the modes until 'revPIP' appears.  
 ⇨ 'Mode' continues to flash.

c. Confirm revPIP mode: Press 'Enter' key.  
 ⇨ The display now shows 'blo' for blow-out.

- d. Preparing for pipetting: Press the pipetting key once to move the piston to its start position.  
 ⇒ The arrow in the display points upwards (aspiration).

- e. Aspirate liquid: Press the pipetting key once to aspirate the liquid.

### NOTICE

When aspirating the liquid, slightly more volume is aspirated than set.

- f. Dispense liquid in revPIP mode: Press the pipetting key once to dispense. On the display, the arrow points down (dispensing). The set volume is now dispensed and some liquid remains in the tip.

- g. Re-aspirate liquid in revPIP mode: Pressing the pipetting key again now resumes the set volume. (Pressing the pipetting key again releases the volume)


- h. Trigger blow-out: After the last pipetting, press 'Enter' key.  
 ⇒ The display now once again shows 'bl' for blow-out.

- i. End pipetting: By pressing the pipetting key once, the blow-out is triggered, and the residual liquid is dispensed.  
 ⇒ After dispensing the residual liquid (blow-out), the display jumps back to the set mode (start position).


## 5.4 Electrophoresis (GEL) mode

Program for loading electrophoresis gels. A previously defined sample volume is aspirated at a high, variable speed and slowly released again.


Volume and speed setting, see Set volume, p. 51 and Set the aspiration and dispensing speed, p. 52.


- a.** Call up menu selection: Pressing the menu key three times takes you to the program selection.  
 ⇒ 'Mode' flashes.


- b.** Set GEL mode: Use one of the arrow keys to scroll through the modes until 'GEL' appears.  
 ⇒ 'Mode' continues to flash.


- c.** Confirm GEL mode: Press 'Enter' key. The display now shows 'blo' for blow-out.

- d.** Preparing for pipetting: Press the pipetting key once to move the piston to its start position.  
 ⇒ The arrow in the display points upwards (aspiration).


- e.** Aspirate liquid: Press the pipetting key once to aspirate the liquid.

- f. In order aspirate more liquid than set (up to max. 110% of the nominal volume), keep the pipetting key pressed during the aspiration process until the desired volume has been aspirated.


⇒ A rhombus appears in the display.

- g. Dispense liquid in GEL mode: To dispense, briefly press the pipetting key once. A rhombus appears in the display. The aspirated volume is slowly released again.


- h. The dispensing of the sample can be interrupted by pressing the pipetting key again.

⇒ The display shows the volume of the liquid dispensed.

- i. Trigger blow-out: After the last pipetting, press 'Enter' key.

⇒ The display now once again shows 'bl o' for blow-out.

- j. End pipetting: By pressing the pipetting key once, the blow-out is triggered, and the residual liquid is dispensed.

⇒ After dispensing the residual liquid (blow-out), the display jumps back to the set mode (start position).


## NOTICE

The GEL mode requires very slow dispensing speeds in order to prevent sample turbulence. In order to ensure optimum dispensing, the dispensing speed is set at the factory. It is considerably slower than the adjustable level 1 and cannot be selected individually.

## 5.5 DISP mode


Program for dispensing an aspirated liquid in partial steps. Slightly more fluid is aspirated than is mathematically necessary.

Volume and speed setting, see Set volume, p. 51 and Set the aspiration and dispensing speed, p. 52.


**a.** Call up menu selection: Pressing the menu key three times takes you to the program selection.

⇒ 'Mode' flashes.


**b.** Set DISP mode: Use one of the arrow keys to scroll through the modes until 'DISP' appears.

⇒ 'Mode' continues to flash.

**c.** Confirm DISP mode: Press 'Enter' key.

⇒ The display now shows 'blo' for blow-out.

**d.** Preparing for pipetting: Press the pipetting key once to move the piston to its start position.

⇒ The arrow in the display points upwards (aspiration).

**e.** Set partial volume: Pressing the arrow key (+/-) sets the volume. Pressing and holding the arrow key causes the volume to change quickly.

⇒ 'VOL' flashes.


**f.** Confirm partial volume: Press 'Enter' key. The display shows the newly set partial volume.

⇒ 'Steps' flashes. The maximum possible number of steps is displayed.


**g.** Set the number of steps: Pressing the arrow key (+/-) sets the number of steps.

⇒ 'Steps' continues to flash.


**h.** Confirm number of steps: Press 'Enter' key.


⇒ The display shows the set number of steps.

**i.** Aspirate liquid: Press the pipetting key once to aspirate the liquid.

**j.** Dispense liquid: Each time the pipetting key is pressed, a dispensing step is performed. The arrow in the display points downwards (dispensing). The 'Step' display shows the number of steps remaining.

**k.** Trigger blow-out: After the last dispensing step, press 'Enter' key.

⇒ The display now once again shows 'blo' for blow-out.

**l.** End dispensing: By pressing the pipetting key once, the blow-out is triggered, and the residual liquid is dispensed.

**m.** After dispensing the residual liquid (blow-out), the display jumps back to the set mode (start position).

## 6 Checking the volume

We recommend testing the device every 3 to 12 months, depending on the level of use. However, the testing cycle can be adapted to meet individual requirements. Gravimetric volume testing of the pipette is carried out according to the following steps and complies with DIN EN ISO 8655, Part 6.

### 1. Setting the nominal volume

- a. Set the maximum specified instrument volume (for procedure, see Pipetting, p. 53).

### 2. Conditioning the pipette

- a. Condition the pipette before testing by aspirating and dispensing the test liquid (distilled water) with a pipette tip five times.

### 3. Performing the test

#### NOTICE

In accordance with DIN EN ISO 8655-2, a tip change is recommended after each individual measurement. An exception to this rule can be made, according to DAkkS guideline DKD-R8-1.

- a. Aspirate the test liquid and pipette into the weighing vessel.
- b. Weigh the pipetted amount with an analysis scale. (refer to the operating manual of the balance manufacturer.)
- c. Calculate the pipetted volume. In doing so, take into account the temperature of the test liquid.
- d. At least 10 pipetting series and weighings in 3 volume ranges (100 %, 50 %, 10 %) are recommended.

### Calculation (for nominal volume)

$x_i$  = weighing results

$n$  = number of weighings

$V_0$  = nominal volume

$Z$  = Correction factor (e.g. 1.0029 µl/mg at 20°C, 1013 hPa)

Mean:

Mean volume:

Accuracy\*:

$$\bar{x} = \frac{\sum x_i}{n}$$

$$\bar{V} = \bar{x} * Z$$

$$A\% = \frac{\bar{V} - V_0}{V_0} * 100$$

Coefficient of variation\*:

Standard deviation\*:

$$CV\% = \frac{100 s}{\bar{V}}$$

$$s = Z * \sqrt{\frac{\sum(x_i - \bar{x})^2}{n - 1}}$$

\* ) = Calculation for accuracy (A%) coefficient of variation (CV%): A% and CV% are calculated using the formulas of statistical quality control.

### NOTICE

Test instructions (SOPs) are available for download at [www.brand.de](http://www.brand.de).

## 7 Accuracy table

Volume range [μl]	Partial volume [μl]	A* ≤ ± %	CV* ≤ %	Sub steps [μl]	Recommended tip type [μl]
0.5-10	10	1.0	0.4	0.01	0.5 - 20
	5	1.5	0.8		
	1	5.0	2.0		
2 - 20	20	1.0	0.4	0.02	0.5 - 20
	10	1.5	0.8		
	2	5.0	2.5		
10 - 200	200	0.8	0.2	0.2	2 - 200
	100	1.2	0.3		
	20	4.0	0.6		
50 - 1000	1000	0.6	0.2	1.0	50 - 1000
	500	1.0	0.3		
	100	3.0	0.6		
250 - 5000	5000	0.6	0.2	5.0	500 - 5000
	2500	1.0	0.3		
	500	3.0	0.6		


\*A = Accuracy, CV = Coefficient of Variation


Final test values based on the nominal volume (= max. volume) printed on the device and the specified partial volumes at the same temperature (20 °C/68 °F) of the device, surroundings and distilled water, in accordance with DIN EN ISO 8655.

# 8 Adjustment – Easy Calibration


## 8.1 Adjustment


The nominal volume or volume to be tested is set, standard mode Pipetting (PIP), e.g. 200.0 µl (for procedure, see PIP mode, p. 55). Ex: Volume according to volume test 201.3 µl.


- a. Bring up CAL mode: Press and hold the MENU key (> 3 sec) until CAL appears.
  - ⇒ The display reads ,off'.
  - ⇒ 'CAL' blinks.


- b. Activate CAL mode: Press one of the arrow buttons to activate the CAL mode.
  - ⇒ The display changes from ,off' to ,on'.
  - ⇒ 'CAL' continues to blink.

- c. Confirm CAL mode: Press the Enter key.
  - ⇒ The display now shows the set pipetting volume.
  - ⇒ 'CAL' blinks.


- d. Setting the volume: Use the arrow keys (+/-) to set the volume, which was previously determined and tested.
  - ⇒ 'CAL' blinks.


- e. Confirm volume: Press the Enter key.
- ⇒ The display shows the tested and corrected volume.
- ⇒ The CAL symbol is continuously displayed to confirm that an adjustment has been made.


## 8.2 Restore factory setting


The CAL symbol permanently shown in the display indicates that an adjustment has been made.

- a. Enter CAL mode: Pressing and holding (> 3 s) the menu key invokes the CAL mode.
- ⇒ The display shows 'on'.
- ⇒ 'CAL' flashes.

- b. Switch off CAL mode: Pressing one of the arrow keys deactivates the CAL mode.
- ⇒ The display changes from 'on' to 'off'.
- ⇒ 'CAL' continues to flash.

- c. Restore factory condition: Press 'Enter' key.
- ⇒ The permanently displayed CAL icon has disappeared.
- ⇒ The instrument has been returned to the original factory condition.

### NOTICE

For the Transferpette® electronic, a volume offset is applied during adjustment (i.e. the volume changes by the same amount over the entire volume range of the pipette). It is recommended to perform the adjustment at 50% of the nominal volume.

### NOTICE

The instrument is permanently adjusted for aqueous solutions but can also be adjusted for solutions of different density, viscosity, and temperature. The Transferpette® electronic can be adjusted in any mode (except GEL mode).

## 9 Disinfection/autoclaving

### 9.1 UV sterilization

The device is resistant to normal exposure to a UV disinfection lamp. The effects of the UV exposure may cause some color change.

### 9.2 Autoclaving

The highlighted part of Transferpette® electronic can be autoclaved at 121°C (250 °F), 2 bar, and with a holding time of at least 15 minutes according to DIN EN 285.

- a. Eject the pipette tip.
- b. Unscrew the pipette shaft from the handle section.
- c. Autoclave the complete pipette shaft without further disassembly.
- d. Allow the pipette shaft to cool down completely and dry.
- e. Screw the pipette shaft back into the handle.
- f. Carry out reference run (rEF).

#### NOTICE

The effectiveness of autoclaving must be verified by the user. Maximum safety is achieved through vacuum sterilization. We recommend the use of sterilization bags.

If the pipette is autoclaved frequently, the piston and seal should be greased with the grease supplied in order to ensure proper movement.

After autoclaving, tighten the connection between the handle and the pipette shaft if necessary.

### 9.3 Reference run (rEF)

After each change of the pipette shaft, a manual reference run must be carried out. The reference run is used for the secure coupling of the piston.

- a. Call up rEF mode: Pressing the 'Menu' and 'Enter' keys simultaneously activates the rEF mode.
  - ⇒ The display shows 'rEF'.
- b. Carry out reference run: Pressing the pipetting key once triggers the reference run.
  - ⇒ A clear functional noise can be heard.
  - ⇒ After the reference run, the display automatically switches back to the previously set program.

# 10 Maintenance


In order to ensure proper functioning, the Transferpette® electronic should be serviced at regular intervals and cleaned as necessary.

## 10.1 Disassembling/cleaning (up to 1,000 µl)

- a. Check the pipette tip cone for damage.
- b. Inspect the piston and seal for contamination.
- c. Check the device for leaks.

We recommend using the leak detector of BRAND, the BRAND PLT unit. As an alternative to this, aspirate a sample and hold the instrument vertically for approx. 10 s. If a drop forms on the pipette tip, refer to Troubleshooting.

### Cleaning

- A Ejector (upper part)
- B Piston unit
- C Seal with spring
- D Shaft
- E Ejector (bottom part)
- S Pipetting shaft

- a. Unscrew the pipette shaft (S) from the handle.
- b. Pull the pipette shaft off the handle. The pipette shaft is attached to the handle by a magnet.
- c. Unscrew the upper part of the ejector unit (A) from the pipette shaft.
- d. Pull out the shaft (B, C and D) from the lower part of the ejector unit (E).
- e. Unscrew the piston unit (B).

#### NOTICE

Do not disassemble the piston unit (B) any further.

- f. Remove the seal with spring (C).
- g. Clean the parts shown with a soap solution or isopropanol, and then rinse with distilled water.
- h. Dry the parts (max. 120 °C/248 °F).
- i. Grease piston and seal with a very thin layer of supplied silicone grease.
- j. Reassemble the cooled parts in reverse order. Only hand-tighten the piston unit and the upper part of the ejector unit (A, B).
- k. Carry out reference run: (rEF), see Reference run (rEF), p. 68.


The individual components shown can be obtained as spare parts, see Order Information, p. 75.

## 10.2 Disassembly/cleaning (250–5000 µl)

- a. Check the pipette tip cone for damage.
- b. Inspect the piston and seal for contamination.
- c. Check the device for leaks.

We recommend using the leak detector of BRAND, the BRAND PLT unit. As an alternative to this, aspirate a sample and hold the instrument vertically for approx. 10 s. If a drop forms on the pipette tip, refer to Troubleshooting.

### Cleaning

- G** Ejector (upper part)
- H** Piston unit
- I** Shaft (bottom part)
- J** Lock
- G'** Ejector (bottom part)
- S** Pipetting shaft

- a. Press side fasteners (J) simultaneously, and pull off lower part of ejector (G').
- b. Unscrew the pipette shaft (H+I) from the handle.
- c. Pull the pipette shaft off the handle. The pipette shaft is attached to the handle by a magnet.
- d. Pull to disconnect the magnetic connection of both components, and remove the upper part of the ejector (G).
- e. Unscrew the piston unit (H) and the lower part of the shaft (I).
- f. Remove the O-ring-seal from the piston unit and clean it.

#### NOTICE

Do not disassemble the piston unit (G) any further.

- g. Clean the piston unit (H) and the bottom part of the shaft (I) with a soap solution or isopropanol, and rinse with distilled water.
- h. Dry the parts (max. 120°C/248°F) and allow them to cool.
- i. Carefully lubricate the inside and outside of the O-ring, and mount it on the piston.
- j. Reassemble the individual components in reverse order.
- k. Carry out reference run: (rEF), see Reference run (rEF), p. 68.

The individual components shown can be obtained as spare parts, see Order Information, p. 75.

## PE filter

PE filter for Transferpette® electronic, 250–5000 µl:

A hydrophobic PE filter is used to prevent liquid from entering the pipette.

Change the filter if it becomes wet or contaminated.

- a. Use a flat object, such as a screwdriver.
- b. Remove the filter carefully, without damaging the tip cone.

Remove the filter before autoclaving!

The device can also be operated without a filter.

## 10.3 Charge and replace rechargeable battery

A fully charged rechargeable battery allows approx. 8 h of continuous pipetting (over 4000 pipetting cycles) of samples of water-like viscosity and density.

### NOTICE

- > Before charging, ensure that the power adapter is suitable for the voltage present in the laboratory.
- > The instrument must not be charged in an explosive environment.
- > The rechargeable battery can be charged only in the Transferpette® electronic.

### Recharge rechargeable battery

- a. Insert the charging cable plug of the power adapter into the socket provided for this purpose at the top of the Transferpette® electronic .
  - ⇒ The charging process starts automatically.
  - ⇒ During the charging process, the bars of the rechargeable battery capacity indicator continually run from the bottom to the top.
  - ⇒ Once the bars of the display have stopped, this means that the rechargeable battery is fully charged.

### Pipetting during the charging process?

During charging, you can continue to work with the Transferpette® electronic . When the rechargeable battery is fully discharged, it takes a few minutes to reach a certain minimum charge capacity. this is necessary for safe operation of the instrument. The last settings made are stored in the EEPROM of the instrument. In the case of complete discharge or when changing the battery, these settings remain saved!


## Replace rechargeable battery

- a.** Open the cover of the battery compartment, remove the rechargeable battery, and pull the plug out of the socket.
- b.** Insert the plug of the new rechargeable battery into the socket and insert the new rechargeable battery.
- c.** Replace the cover of the battery compartment and close it.

Remove the rechargeable battery from the instrument during longer breaks in operation.

## Battery indicator after reinserting a battery pack

After inserting a rechargeable battery, the full capacity indicator appears in the display with a flashing frame (the instrument does not yet recognize the charge status at first). After 3.5 h charging time – safe complete charging of the battery – the frame stops flashing.

### NOTICE

Always charge for at least 3.5 h after inserting a battery! The full charge capacity is reached after several charge/discharge cycles!



## 10.4 Rechargeable battery regeneration function

### Refresh function

To extend the service life and increase the performance of the batteries, the Transferpette® electronic has a regeneration function (refresh function). This function allows the batteries to be fully discharged and recharged under program control. To optimize the performance of the batteries, the refresh function should be used from time to time.

## Carry out refresh function

- a. Insert the charging cable plug (connection) of the power adapter into the socket provided for this purpose at the top of the Transferpette® electronic .
- b. Press the lower arrow key for longer than 3 s. During discharging, the capacity bars of the battery display run continuously from top to bottom.
- c. After discharging (up to 3 h), the charging process (3.5 h) is started automatically. During charging, the capacity bars of the battery display run continuously from bottom to top.

## Canceling the refresh function


Pressing any key terminates the program. The instrument automatically switches back to the standard pipetting mode (PIP) and to the nominal volume, and the normal loading process is automatically started, see Charge and replace rechargeable battery, p. 71. Unplugging the power supply also terminates the program. Cancellation of the refresh function must not be performed at the end of the discharge cycle.

# 11 Troubleshooting

Fault	Display	Cause	Corrective action
Instrument does not respond	ERR 1	Rechargeable battery empty or defective	Charge the battery for at least 5 min without actuation. Then continue working with the charging cable only until the battery is recharged. Replace the rechargeable battery if necessary
		Electronic components defective	Send instrument for repair
Instrument does not respond	ERR 2	Electronic components defective	Send instrument for repair

Fault	Display	Cause	Corrective action
Instrument does not respond	ERR 3	Unforeseen program error	Error confirmation by pressing the 'Enter' key, instrument is reinitialized
Instrument does not respond	ERR 4	No rechargeable battery in the instrument	Insert rechargeable battery
		Rechargeable battery defective	Replace rechargeable battery
		Electronic components defective	Send instrument for repair
Tip dripping/instrument leaking or volume error	—	Unsuitable tip	Use only high-quality tips
		Tip not seated tightly	Press tip on more firmly/other interchangeable clip
		Piston, shaft, or seal dirty or damaged	Clean instrument/replace seal, grease piston
No indication in the display	—	Electrostatic discharge	Remove rechargeable battery and insert again
		Electronic components defective	Send instrument for repair
Aspiration is not possible	—	Motor has no connection to the pipetting unit.	Carry out reference run (rEF), see Reference run (rEF), p. 68.

## 12 Product markings

Symbol or number	Meaning
	With this mark, we confirm that the product complies with the requirements set out in the EC Directives and has been subjected to the specified testing procedures.
	UKCA: United Kingdom Conformity Assessed With this mark, we confirm that the product meets the requirements specified in the UK Designated Standards .
	The instrument is marked in accordance with the German Measurement and Calibration Act as well as the Measurement and Calibration Regulation. Character sequence DE-M (DE for Germany), framed by a rectangle, as well as the two last digits of the year the marking was affixed.
www.brand.de/ip	Patent information
XXZXXXXX	Serial number
	Follow the instructions listed on the instrument, the accessory parts and in the operating manual.

Symbol or number	Meaning
	The instrument or rechargeable battery should be disposed of properly.
	China RoHS (EFUP) EFUP defines the time period in years, in which the hazardous materials found in the electrical and electronic equipment should not leak or mutate under normal operating conditions. When used under normal conditions, such electrical or electronic products do not lead to severe environmental pollution, severe personal injuries or damage to the user's property.
	This electrical instrument may not be disposed of in household waste.

## 13 Order Information

### 13.1 Ordering Information

#### Transferpette® electronic

Volume	0.5–10 µl	2–20 µl	10–200 µl	50–1000 µl	250–5000 µl
with power adapter AC 100–240 V ~50–60 Hz	Order no.	Order no.	Order no.	Order no.	Order no.
for Europe	705299	705300	705303	705306	705307
for the UK/Ireland	705309	705310	705313	705316	705317
for the USA/ Japan	705319	705320	705323	705326	705327
for Australia	705329	705330	705333	705336	705337
Without power adapter	705339	705340	705343	705346	705347

#### Power adapters (AC 100–240 V~50–60 Hz)

Description	Order no.
for Europe	705350
for the UK/Ireland	705351
for the USA/Japan	705352
for Australia	705353

## Spare battery

Description	Order no.
Spare battery for Transferpette® electronic	705500

## Silicone grease

Description	Order no.
Silicone grease for Transferpette® electronic up to 1,000 µl	705502
Silicone grease for Transferpette® electronic 250–5000 µl	703677

## PLT unit

Description	Order no.
PLT unit (pipette leak testing unit)	703970

## Charging stand with power adapter

Description	Order no.
Charging stand with power adapter (AC 100–240 V ~50–60 Hz) for 3 Transferpette® electronic up to 1000 µl.	
for Europe	705390
for the UK/Ireland	705391
for the USA/Japan	705392
for Australia	705393

## Quality pipette tips from BRAND

Volume	Material	Packaged unit	Order no.
		Packed in bulk	
0.1–20 µl	PP	2000	732002
0.5–20 µl	PP	2000	732004
1–50 µl	PP	2000	732006
2–200 µl	PP	1000	732008
50–1000 µl	PP	1000	732012
5 ml	PP	200	702595


## Filter, 5 ml

Description	Order no.
Filter for Transferpette® electronic 5 ml, packing unit 25 pieces	704652

## 13.2 Spares

### 13.2.1 Transferpette® electronic to 1000 µl

The design and dimensions of the spare parts correspond to the respective nominal volume. (Fig. Spare parts for Transferpette® electronic 50–1000 µl.)

A Ejector (upper part)

B Piston unit

C Spring with seal


D Shaft with ejector spring

E Ejector (bottom part)

S Complete shaft

Volume	A	B	C	D	E	L
0.5–10 µl	705510	705518	–	705538	705548	705528
2–20 µl	705510	705520	705530	705539	705550	705529
10–200 µl	705510	705523	704662	705537	705553	705546
50–1000 µl	705510	705525	704663	705541	705555	705547

## 13.2.2 Transferpette® electronic 250–5000 µl

- G** Ejector (upper part)
- H** Piston unit
- S** Complete shaft

- G'** Ejector (bottom part)
- I** Shaft (bottom part)

Volume	G + G'	H	I	S
250 - 5000 µl	7299	705526	703247	705549

## 13.3 Additional accessories

### 13.4 Accessories for Transferpette® electronic 10–200 µl and 50–1000 µl

#### Interchangeable clips

Attachable interchangeable clips (spacers) and pipette shafts with flexible tip mounting cone enable optimum fitting accuracy and lowest possible ejection forces with most commercially available pipette tips.

Volume	Packaged unit	Order No.
10–200 µl	3	705560
50–1000 µl	3	705562

## Pipette shaft soft (complete)

Pipette shafts complete for Transferpette® electronic , single channel, soft shaft

Volume	Order No.
10–200 µl	705543
50–1000 µl	705545

## 14 Repairs

### 14.1 Sending for repair

#### NOTICE

Transporting of hazardous materials without a permit is a violation of federal law.

#### Clean the instrument thoroughly and decontaminate!

- When returning products, please enclose a general description of the type of malfunction and the media used. If information regarding media used is missing, the instrument cannot be repaired.
- Only send the device without a battery installed.
- Shipment is at the risk and the cost of the sender.

#### Outside USA and Canada

Complete the “Declaration on Absence of Health Hazards” and send the instrument to the manufacturer or supplier. Ask your supplier or manufacturer for the form. The form can also be downloaded from [www.brand.de](http://www.brand.de).

#### Within USA and Canada

Please clarify the requirements for the return delivery with BrandTech Scientific, Inc **before** sending the instrument in for service.

Return only cleaned and decontaminated instruments to the address provided with the Return Authorization Number. Place the Return Authorization number so that it is clearly visible on the outside of the package.

## Contact addresses

### Germany:

BRAND GMBH + CO KG  
Otto-Schott-Str. 25  
97877 Wertheim (Germany)  
T +49 9342 808 0  
F +49 9342 808 98000  
info@brand.de  
www.brand.de

### USA and Canada:

BrandTech<sup>®</sup> Scientific, Inc.  
11 Bokum Road  
Essex, CT 06426-1506 (USA)  
T +1-860-767 2562  
F +1-860-767 2563  
info@brandtech.com  
www.brandtech.com

### India:

BRAND Scientific Equipment Pvt. Ltd.  
303, 3rd Floor, 'C' Wing, Delphi  
Hiranandani Business Park,  
Powai  
Mumbai-400 076 (India)  
T +91 22 42957790  
F +91 22 42957791  
info@brand.co.in  
www.brand.co.in

### China:

BRAND (Shanghai) Trading Co., Ltd.  
Guangqi Culture Plaza  
Room 506, Building B  
No. 2899, Xietu Road  
Shanghai 200030 (P.R. China)  
T +86 21 6422 2318  
F +86 21 6422 2268  
info@brand.com.cn  
www.brand.cn.com

## 15 Calibration service

The ISO 9001 and GLP guidelines require regular inspection of your volume measuring devices. We recommend performing a volume check every 3 to 12 months. The cycle is dependent on the individual requirements of the device. Checks should be performed more frequently, in case of high frequency of use or the use of aggressive media.

The detailed testing instructions are available for download on [www.brand.de](http://www.brand.de) or [www.brandtech.com](http://www.brandtech.com).

BRAND also offers the possibility to have your devices calibrated by our factory calibration service or by the BRAND DAkkS laboratory. Simply send us the device to be calibrated, accompanied by details about which type of calibration you would like. The device will be returned to you after a few days together with a test report (factory calibration) or a DAkkS calibration certificate. More information can be obtained from your dealer or directly from BRAND. The order document is available for download on [www.brand.de](http://www.brand.de) (see Technical Documents).

## 16 Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operating or unauthorized repairs of the device or for the consequences of normal wear and tear, especially of wearing parts such as pistons, seals, valves and the breakage of glass. The same applies for failure to follow the instructions of the operating manual. We are not liable for damage resulting from disassembly beyond that described in the operating manual or if non-original spare parts or components have been installed.

### USA and Canada:

Find more warranty information on [www.brandtech.com](http://www.brandtech.com).

## 17 Disposal

This symbol means that at the end of their service life, batteries/accumulators and electronic devices must be disposed of separately from household waste (unsorted municipal waste).

Electronic devices must be disposed of in accordance with Directive 2012/19/ EU of the European Parliament and of the Council from July 04, 2012 on waste from electrical and electronic equipment and in compliance with national disposal regulations.

Both batteries and accumulators (rechargeable batteries) contain materials that can be damaging to the environment and human health. Therefore, they must be properly disposed of in accordance with Directive 2006/66/EC of the European Parliament and of the Council from September 06, 2006 on batteries and accumulators and in compliance with national disposal regulations. Only dispose of fully discharged batteries and accumulators.

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